


PATCHES

2-3 COMPARTMENTS

Regulatory category Medical devices/others

Indication Pain relief

Unit size 2 patches of 2-3 compartments



geiserpharma



ACTIVE INGREDIENT (S)

Iron powder
Activated carbon

UNIT SIZE (DETAILS)

2 patches of 22 x 10 cm display
in 3 compartments
2 patches of 14 x 10 cm display
in 2 compartments

PRIMARY PACKAGING

Plastic bag

SECONDARY PACKAGING

Folding carton box

DOSSIER

Technical/ quality pharma
dossier, English language.
Cytotoxicity test (adhesive
layer), irritation test and sensi-
bility test included

PRODUCT DESCRIPTION

Warming patch to provide pain relief by propor-
tioning 8 hours of constant and penetrating
warmth when placed on contracted muscles.
The warmth generated by activated carbon &
iron powder natural thermal ingredients, stimu-
lates blood circulation, allowing a blood flow
enriched with oxygen and nutrients which
relaxes tense and contracted muscles. The
patch needs to be applied directly on the skin
using for this purpose an adhesive appropriate
for using on human skin which assures perfect
adhesion of the patch during min. 8 hours. Easy
application, no adhesive remains. Used on lower
back area, shoulders or neck area.

TECHNOLOGY FOR HEALTH

PHARMACEUTICAL DEVELOPMENT

T. +34 948 101 943

geiserpharma@geiserpharma.com

www.geiserpharma.com